

Distributors Info:

SpotCheck™

Rapid Food Residue Test

Designed for companies looking for an accurate, easy-to-use, non-instrument based hygiene surface test, SpotCheck offers an inexpensive way to monitor the effectiveness of cleaning programs.

SpotCheck can detect the presence of glucose in 60 seconds or less. Glucose is present in 85% of food items. Left on contact surfaces, food residues provide nutrients for bacterial growth which could lead to dangerous contamination and a diminished shelf life. A quick color change from clear to green of the self-contained reagent indicates that glucose residue is present.




IF IT TURNS GREEN, IT'S NOT CLEAN!

Key Features and Benefits	
✓	Results in under 60 seconds
✓	Sensitive to low levels of contamination
✓	Easy to interpret
✓	Semiquantitative results
✓	No instrumentation
Ideal for	
✓	Bakeries
✓	Grocery stores
✓	Fast-food stores and restaurants
✓	Fruit and vegetable processors
✓	Beverage manufacturers
✓	Dressing and sauce processors
✓	Confectioneries

NO INSTRUMENTATION

RESULTS IN 60 SEC



Americas
 941 Avenida Acaso
 Camarillo, CA 93012
 USA
Tel: 805.388.8007
Fax: 805.388.5531
Email: enquiries@hygiene.net
Web: www.hygiene.net

Europe
 Unit 25 Wenta Business Centre
 Colne Way Watford
 Hertfordshire WD2 4ND UK
Tel: +44 (0)1923 818821
Fax: +44 (0)1923 81825
Email: enquiries@hygiene.net
Web: www.hygiene.net

ORDER # L0100

(100 tests per box)

Distributors Info:

SpotCheck Plus™

Rapid Food Residue Test with Added Detection

SpotCheck Plus has the same detection capability of SpotCheck with the additional benefit of lactose detection. Simply swab a surface, release the reagent and if any residue containing glucose and/or lactose is present the reagent will turn green, providing a semiquantitative measure of the surface hygiene. The more contamination present, the quicker the color change and darker the color. SpotCheck Plus instantly validates the hygiene of a surface, allowing immediate corrective action to be taken if necessary.




IF IT TURNS GREEN, IT'S NOT CLEAN!

Key Features and Benefits	
✓	Results in under 60 seconds
✓	Sensitive to low levels of contamination
✓	Easy to interpret
✓	Semiquantitative results
✓	No instrumentation
Ideal for	
✓	Dairies and creameries
✓	Bakeries
✓	Grocery stores
✓	Fast-food stores and restaurants
✓	Fruit and vegetable processors
✓	Beverage manufacturers
✓	Dressing and sauce processors
✓	Confectioneries

NO INSTRUMENTATION

RESULTS IN 60^{SEC}



Americas
 941 Avenida Acaso
 Camarillo, CA 93012
 USA
Tel: 805.388.8007
Fax: 805.388.5531
Email: enquiries@hygiena.net
Web: www.hygiena.net

Europe
 Unit 25 Wenta Business Centre
 Colne Way Watford
 Hertfordshire WD2 4ND UK
Tel: +44 (0)1923 818821
Fax: +44 (0)1923 81825
Email: enquiries@hygiena.net
Web: www.hygiena.net

ORDER # GL100

(100 tests per box)